


# Blood Coagulator



**Function:** To coagulate continuously the raw blood flow with injected direct steam into a blood clot which is ready for further processing/drying.

**Features:** Full stainless steel execution with direct steam injection points and sight glasses in the coagulator discharge chute. The stainless steel coagulator is placed on a support frame.

**Capacity:** 1- 7 ton of raw blood per hour

**MAVITEC**

Galileistraat 32, 1704 SE Heerhugowaard, The Netherlands  
Tel.: +31-72-574 59 88 - Fax: +31-72-574 55 48  
E-mail: [info@mavitec.com](mailto:info@mavitec.com) - [www.mavitec.com](http://www.mavitec.com)

**DUPPS**

Germantown, Ohio U.S.A.  
Tel: +1-937-855 6555 • Fax: +1-937-855 6554  
E-mail: [info@dupps.com](mailto:info@dupps.com) • [www.dupps.com](http://www.dupps.com)